

Human RNASE1/RIB1/RNS1 (His tag) recombinant protein



Catalog Number: 518085

General Information

Gene Name Synonym

Ribonuclease Pancreatic, HP-Rnase, RIB-1, RNase UpI-1, Ribonuclease 1, RNase 1, Ribonuclease A, RNase A, RNASE1, RIB1, RNS1

Protein Construction

A DNA sequence encoding the Human RNASE1/RIB1/RNS1 Lys29-Thr156 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 10% Glycerol, pH 7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human RNASE1/RIB1/RNS1 (His tag)
recombinant protein