

SNX7 antibody

Catalog Number: 115485

Product name

SNX7 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SNX7 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 45kd, 52kd

Preparation

This antibody was obtained by immunization of SNX7 recombinant protein (Accession Number: XM_017001428). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115485(SNX7 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:115485(SNX7 antibody) at dilution of 1:50 (under 10x lens)