

# COX16 antibody



Catalog Number: 109485

---

## Product name

COX16 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

COX16 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, A375 cells, HepG2 cells, human heart tissue. Observed molecular weight by Western-blot: 12kd

## Preparation

This antibody was obtained by immunization of COX16 recombinant protein (Accession Number: NM\_016468). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:109485(COX16 antibody) at dilution of 1:1000