

# GATAD2B antibody



Catalog Number: 110885

## Product name

GATAD2B antibody

## Immunogen

[Human GATAD2B Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

GATAD2B Rabbit Polyclonal antibody. Positive IF detected in A549 cells. Positive WB detected in A549 cells, HEK-293 cells, HeLa cells, K-562 cells, MCF-7 cells, PC-3 cells. Observed molecular weight by Western-blot: 69 kDa

## Preparation

This antibody was obtained by immunization of GATAD2B recombinant protein (Accession Number: XM\_005245364). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

## Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:110885(GATAD2B Antibody) at dilution of



Immunofluorescent analysis of A549 cells using Catalog No:110885(GATAD2B Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)