

JIP3 antibody



Catalog Number: 111885

Product name

Polyclonal

JIP3 antibody

Ig Type

Specificity

Rabbit IgG

Human, Rat, Mouse; other species not tested.

Applications

Antibody description

ELISA, WB

JIP3 Rabbit Polyclonal antibody. Positive WB detected in rat brain tissue. Observed molecular weight by Western-blot: 147 kDa

Dilutions

Recommended Dilution:


Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of JIP3 recombinant protein (Accession Number: NM_015133). Purification method: Antigen affinity purified.

Validations

Formulation



PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111885(JIP3 Antibody) at dilution of 1:600

Store at -20°C. DO NOT ALIQUOT

Clonality