

Anti-phospho-Rb (Ser780) Rabbit antibody



Catalog Number: 167885

Product name

Anti-phospho-Rb (Ser780) Rabbit antibody

Specificity

Human

Antibody description

phospho-Rb (Ser780) Rabbit monoclonal antibody

Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Ser780 of human Rb

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions

WB: 1/2000

IHC: 1/20

Validations



Western blot detection of Rb (phospho-Ser780) in HeLa, A549, HL-60, U2OS, C6 cell lysates using Rb (phospho-Ser780) Rabbit mAb(1:1000 diluted). Predicted band size: 106kDa. Observed band size: 110kDa.