

GSTM1 antibody



Catalog Number: 111185

Product name

GSTM1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

GSTM1 Rabbit Polyclonal antibody. Positive WB detected in mouse ovary tissue, HeLa cells, mouse kidney tissue, mouse liver tissue. Positive IHC detected in human gliomas tissue. Observed molecular weight by Western-blot: 27 kDa

Preparation

This antibody was obtained by immunization of GSTM1 recombinant protein (Accession Number: NM_146421). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



mouse ovary tissue were subjected to SDS PAGE followed by western blot with Catalog No:111185(GSTM1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:111185(GSTM1 antibody) at dilution of 1:50 (under 10x lens)