

NINJ2 antibody



Catalog Number: 113185

Product name

Polyclonal

NINJ2 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

NINJ2 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

Dilutions

Preparation

Recommended Dilution:


This antibody was obtained by immunization of NINJ2 recombinant protein (Accession Number: BC057766). Purification method: Antigen affinity purified.

WB: 1:500-1:5000

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Validations



Storage

Store at -20°C. DO NOT ALIQUOT

Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113185(NINJ2 Antibody) at dilution of 1:1000

Clonality