

PREPL antibody

Catalog Number: 114185

Product name

PREPL antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PREPL Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 72 kDa

Preparation

This antibody was obtained by immunization of PREPL recombinant protein (Accession Number: BC013193). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114185(PREPL antibody) at dilution of 1:400