

# PPIC Polyclonal Antibody



Catalog Number: 166185

## Product name

PPIC Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to PPIC

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 23-212 of human PPIC (NP\_000934.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:100

## Validations



Western blot - PPIC Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using PPIC antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.



Immunohistochemistry - PPIC Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human prostate using PPIC antibody at dilution of 1:100 (40x lens).