

# Anti-MIF Rabbit antibody



Catalog Number: 167185

## Product name

Anti-MIF Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

MIF Rabbit polyclonal antibody

## Preparation

Antigen: Recombinant protein of human MIF

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, FC


## Dilutions

WB: 1/1000

ICC/IF: 1/20

FC: 1/100

## Validations



Western blot detection of MIF in C6,3T3 cell lysates using MIF Rabbit pAb(1:1000 diluted).Predicted band size:13kDa.Observed band size:13kDa.