

# CXCL8/IL8 antibody



Catalog Number: 107185

## Product name

CXCL8/IL8 antibody

## Specificity

Human; other species not tested.

## Antibody description

CXCL8/IL8 Mouse Monoclonal antibody. Positive WB detected in ag10552. Observed molecular weight by Western-blot: 37 kDa

## Preparation

This antibody was obtained by immunization of CXCL8/IL8 recombinant protein (Accession Number: NM\_000584). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:2000-1:20000

## Validations



Recombinant protein were subjected to SDS PAGE followed by western blot with Catalog No:107185 at dilution of 1:1000