

Ozanezumab reference antibody



Catalog Number: 177185

Product name

Ozanezumab reference antibody

Specificity

Human RTN4 / NOGO

Antibody description

Anti-RTN4 / NOGO Antibody (Ozanezumab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG1

Applications


ELISA

Validations



SDS-PAGE

Anti-RTN4 / NOGO Antibody (ozanezumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.



SEC-HPLC

The purity of Anti-RTN4 / NOGO Antibody (ozanezumab) is 96.23%, determined by SEC-HPLC.