

ICA1L antibody



Catalog Number: 111585

Product name

ICA1L antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

ICA1L Rabbit Polyclonal antibody. Positive IF detected in MCF-7.

Preparation

This antibody was obtained by immunization of ICA1L recombinant protein (Accession Number: NM_001288624). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IF

Dilutions

Recommended Dilution:

IF: 1:10-1:100

Validations



Immunofluorescent analysis of MCF-7 cells, using ICA1L antibody Catalog No:111585 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).