

Mouse Fas Ligand / FASLG / CD95L (His Tag) recombinant protein



Catalog Number: 502585

General Information

Protein Construction

A DNA sequence encoding the mouse Faslg (AAB02915.1) (Pro132-Leu279) was expressed with a polyhistidine tag at the N-terminus.

Organism

Mouse

Expression Host

Yeast

QC Testing

Purity

> 90 % as determined by SDS-PAGE.

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant mouse Faslg consists of 158 amino acids and predicts a molecular mass of 18.2

kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse Fas Ligand / FASLG / CD95L Protein (His Tag)