

# Anti-non-muscle Myosin IIA Rabbit antibody



Catalog Number: 167585

## Product name

Rabbit IgG

Anti-non-muscle Myosin IIA Rabbit antibody

## Applications

## Specificity

WB, IHC-P, IP

Human, Mouse, Rat

## Dilutions

## Antibody description

WB: 1/1000

non-muscle Myosin IIA Rabbit monoclonal antibody

IHC: 1/50

IP: 1/20


## Preparation

Antigen: A synthetic peptide of human non-muscle Myosin IIA

## Validations

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA



## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

Western blot detection of non-muscle Myosin IIA in HeLa, A549, HL-60, U2OS, C6 cell lysates using non-muscle Myosin IIA Rabbit mAb(1:1000 diluted). Predicted band size: 227kDa. Observed band size: 227kDa.

## Ig Type