

SFRP2 antibody



Catalog Number: 107585

Product name

SFRP2 antibody

Specificity

Human, Rat; other species not tested.

Antibody description

SFRP2 Mouse Monoclonal antibody. Positive WB detected in rat brain tissue, COLO 320 cells, Y79 cells. Observed molecular weight by Western-blot: 38 kDa

Preparation

This antibody was obtained by immunization of SFRP2 recombinant protein (Accession Number: NM_003013). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107585(SFRP2 Antibody) at dilution of 1:2000