

GJB6 antibody

Catalog Number: 110985

Product name

GJB6 antibody

Immunogen

[Human GJB6 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

GJB6 Rabbit Polyclonal antibody. Positive IHC detected in human oesophagus tissue, human skin cancer tissue.

Preparation

This antibody was obtained by immunization of GJB6 recombinant protein (Accession Number: NM_001110221). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC


Dilutions

Recommended Dilution:
IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human oesophagus tissue slide using Catalog No:110985(GJB6 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human oesophagus tissue slide using Catalog No:110985(GJB6 Antibody) at dilution of 1:50 (under 40x lens)