

Anti-LAMP1 Rabbit antibody



Catalog Number: 166985

Product name

Anti-LAMP1 Rabbit antibody

Specificity

Human

Antibody description

LAMP1 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human LAMP1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions

WB: 1/1000

IHC: 1/50

Validations



Western blot detection of LAMP1 in HeLa cell lysates using LAMP1 Rabbit pAb(1:1000 diluted). Predicted band size:45kDa.Observed band size:90-120kDa.