

XIAP antibody

Catalog Number: 116985

Product name

XIAP antibody

Specificity

Human, Rat; other species not tested.

Antibody description

XIAP Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, C6 cells, HEK-293 cells, HeLa cells, Tunicamycin treated HepG2 cells. Observed molecular weight by Western-blot: 45-50 kDa

Preparation

This antibody was obtained by immunization of XIAP recombinant protein (Accession Number: NM_001204401). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116985(XIAP



antibody) at dilution of 1:1000

Tunicamycin treated HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:116985(XIAP Antibody) at dilution of 1:600