

# XIAP antibody

Catalog Number: 116985

## Product name

XIAP antibody

## Specificity

Human, Rat; other species not tested.

## Antibody description

XIAP Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, C6 cells, HEK-293 cells, HeLa cells, Tunicamycin treated HepG2 cells. Observed molecular weight by Western-blot: 45-50 kDa

## Preparation

This antibody was obtained by immunization of XIAP recombinant protein (Accession Number: NM\_001204401). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116985(XIAP



antibody) at dilution of 1:1000 Tunicamycin treated - +

Tunicamycin treated HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:116985(XIAP Antibody) at dilution of 1:600