

# Anti-IRAK4 Rabbit antibody



Catalog Number: 168985

## Product name

Anti-IRAK4 Rabbit antibody

## Specificity

Human

## Antibody description

IRAK4 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human IRAK4

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, FC

## Dilutions


WB: 1/1000-1/5000

IHC: 1/20-1/500

ICC/IF: 1/20

FC: 1/20

## Validations



Western blot detection of IRAK4 in Jurkat, HeLa cell lysates using IRAK4 Rabbit mAb(1:1000 diluted). Predicted band size: 52kDa. Observed band size: 52kDa.