

# Anti-NFI-C/CTF (4E10) Mouse antibody



Catalog Number: 168285

## Product name

Anti-NFI-C/CTF (4E10) Mouse antibody

## Specificity

Human

## Antibody description

NFI-C/CTF (4E10) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human NFIC protein fragments expressed in E.coli.

## Formulation

Ascites containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications


WB, ICC/IF

## Dilutions


WB: 1/200

;ICC: 1/200

## Validations



Immunocytochemistry staining of HeLa cells fixed with 4% Paraformaldehyde and using anti-NFIC



mouse mAb (dilution 1:200).

Western blot detection of NFIC in HeLa NE cell lysates using NFIC mouse mAb (1:200 diluted). Predicted band size: 56kDa. Observed band size: 56kDa.