

Mouse IL36G (His Tag) recombinant protein



Catalog Number: 502685

General Information

Gene Name Synonym

Interleukin-1 family member 9

Protein Construction

A DNA sequence encoding the mouse IL36g (Q8R460) (Gly13-Ser164) was expressed with a polyhistidine tag at the N-terminus.

Organism

Mouse

Expression Host

E. coli

QC Testing

Purity

> 85 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant mouse IL36g consists of 168

amino acids and predicts a molecular mass of 19.4 KDa. It migrates as an approximately 20 KDa band in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse IL36G / IL1F9 Protein (His Tag) SDS-PAGE