

GNA15 Polyclonal Antibody



Catalog Number: 163685

Product name

GNA15 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to GNA15

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-374 of human GNA15 (NP_002059.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations




Western blot - GNA15 Polyclonal Antibody

Western blot analysis of extracts of mouse heart, using GNA15 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.



Immunohistochemistry - GNA15 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human tonsil using GNA15 antibody at dilution of 1:100 (40x lens).



GNA15 Polyclonal Antibody



Catalog Number: 163685

Immunohistochemistry - GNA15 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human breast cancer using GNA15 antibody at dilution of 1:100 (40x lens).