

Anti-Aurora B Rabbit antibody



Catalog Number: 166685

Product name

Anti-Aurora B Rabbit antibody

Specificity

Human

Antibody description

Aurora B Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human Aurora B

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, IP

Dilutions


WB: 1/1000-1/2000

IHC: 1/50-1/200


ICC/IF: 1/50-1/200

IP: 1/50


Validations



Immunohistochemical analysis of paraffin-embedded human tonsil, using Aurora B Antibody.



Western blot analysis of Aurora B expression in Jurkat and HeLa cell lysates using Aurora B antibody at 1/1000 dilution. Predicted band size: 39KDa. Observed band size: 39KDa.



Western blot analysis of Aurora B expression in HeLa cell lysate.