

# Anti-phospho-JNK1 (Thr183/Tyr185) Rabbit antibody



Catalog Number: 167685

## Product name

Anti-phospho-JNK1 (Thr183/Tyr185) Rabbit antibody

## Specificity

Human

## Antibody description

phospho-JNK1 (Thr183/Tyr185) Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Thr183/Tyr185 of human JNK1

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, IP

## Dilutions

WB: 1/1000

IP: 1/20

## Validations



Western blot detection of phospho-JNK1 (Thr183/Tyr185) in HeLa, A549, C6 cell lysates using phospho-JNK1 (Thr183/Tyr185) Rabbit mAb(1:1000 diluted). Predicted band size: 48,53kDa. Observed band size: 46,54kDa.