

# Anti-OXGR1/GPR99 Rabbit antibody



Catalog Number: 168685

## Product name

Anti-OXGR1/GPR99 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

OXGR1/GPR99 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human OXGR1

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 1/1000-1/5000

## Validations



Western blot detection of OXGR1/GPR99 in HeLa cell lysates using OXGR1/GPR99 Rabbit mAb(1:1000 diluted). Predicted band size:38kDa. Observed band size:38kDa.