

IL-17F antibody



Catalog Number: 107385

Product name

Monoclonal

IL-17F antibody

Ig Type

Specificity

Mouse IgG1

Human, Mouse; other species not tested.

Applications

Antibody description

ELISA, WB

IL-17F Mouse Monoclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 37 kDa

Dilutions


Recommended Dilution:

Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of IL-17F recombinant protein (Accession Number: NM_052872). Purification method: Protein G purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107385(IL17F antibody) at dilution of 1:1000

Clonality