

CLEC4G antibody

Catalog Number: 109385

Product name

CLEC4G antibody

Specificity

Human; other species not tested.

Antibody description

CLEC4G Rabbit Polyclonal antibody. Positive WB detected in HUVEC cells, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 40 kDa

Preparation

This antibody was obtained by immunization of CLEC4G recombinant protein (Accession Number: NM_198492). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Validations



HUVEC cells were subjected to SDS PAGE followed by western blot with Catalog No:109385(CLEC4G Antibody) at dilution of 1:300



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:109385 (CLEC4G Antibody) at dilution of 1:1000