

NLRP7 Polyclonal Antibody



Catalog Number: 160785

Product name

NLRP7 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to NLRP7

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 800-1000 of human NLRP7 (NP_631915.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:100

Validations



Western blot - NLRP7 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using NLRP7 Antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 40s.