

Human RAB31/RAB22B (His Tag) recombinant protein



Catalog Number: 501785

General Information

Gene Name Synonym

Ras-related protein Rab-22B

Protein Construction

A DNA sequence encoding the human RAB31 (NP_006859.2) (Met1-Cys195) was expressed with a polyhistidine tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant human RAB31 consists of 213

amino acids and predicts a molecular mass of 23.9 kDa. It migrates as an approximately 26 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 50 mM Tris, pH 8.0.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C.

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human RAB31 / Rab22B Protein (His Tag) SDS-PAGE