

# PGK2 antibody

Catalog Number: 113785

---

## Product name

PGK2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PGK2 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, HeLa cells, HepG2 cells, human heart tissue, human testis tissue, mouse brain tissue. Observed molecular weight by Western-blot: 45 kDa

## Preparation

This antibody was obtained by immunization of PGK2 recombinant protein (Accession Number: NM\_138733). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:113785(PGK2 antibody) at dilution of 1:1200