

RON, MST1R antibody



Catalog Number: 114785

Product name

RON, MST1R antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RON, MST1R Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue. Observed molecular weight by Western-blot: 140 kDa

Preparation

This antibody was obtained by immunization of RON, MST1R recombinant protein (Accession Number: BQ936677 NM_002447). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:114785(RON, MST1R antibody) at dilution of 1:300