

# ATXN1 Polyclonal Antibody



Catalog Number: 160095

## Product name

ATXN1 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to ATXN1

## Preparation

Antigen: A synthetic peptide of human ATXN1

## Formulation

PBS with 0.02% sodium azide, pH7.3.

## Storage

Store at 4°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:200 - 1:500

IF 1:20 - 1:50

## Validations



Western blot - ATXN1 Polyclonal Antibody

Western blot analysis of extracts of CEM cells, using ATXN1 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.