

# HOXB7 Polyclonal Antibody



Catalog Number: 165095

## Product name

HOXB7 Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to HOXB7

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human HOXB7 (NP\_004493.3).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:100

## Validations



Western blot - HOXB7 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using HOXB7 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure



Immunofluorescence - HOXB7 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using HOXB7 antibody.