

A1BG antibody

Catalog Number: 107095

Product name

A1BG antibody

Specificity

Human; other species not tested.

Antibody description

A1BG Mouse Monoclonal antibody. Positive WB detected in human plasma tissue, human blood tissue, human heart tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 72 kDa

Preparation

This antibody was obtained by immunization of A1BG recombinant protein (Accession Number: BC035719). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:2000-1:20000

IF: 1:50-1:500

Validations



human plasma tissue were subjected to SDS PAGE followed by western blot with Catalog No:107555(A1BG Antibody) at dilution of 1:8000



Immunofluorescent analysis of (-20oc Ethanol) fixed HeLa cells using Catalog No:107555(A1BG Antibody) at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)