

Anti-SOCS1 Rabbit antibody



Catalog Number: 169095

Product name

Anti-SOCS1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

SOCS1 Rabbit polyclonal antibody

Preparation

Antigen: The antiserum was produced against synthesized peptide derived from human SOCS-1. AA range:49-98

Formulation

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ICC/IF, WB, IHC-P, ELISA


Dilutions

IF: 1/50-200

WB: 1/500-2000; ELISA 1/10000-20000

IHC: 1/50-300

Validations



Western Blot analysis of Jurkat cell lysate using SOCS-1 Polyclonal Antibody