

# INHA Polyclonal Antibody



Catalog Number: 163495

## Product name

INHA Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to INHA

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 19-366 of human INHA (NP\_002182.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

## Validations



Western blot - INHA Polyclonal Antibody

Western blot analysis of extracts of rat testis, using INHA antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunohistochemistry - INHA Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using INHA antibody at dilution of 1:200 (40x lens).