

# RAC2 antibody

Catalog Number: 107495

## Product name

RAC2 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

RAC2 Mouse Monoclonal antibody. Positive IHC detected in human spleen tissue, human tonsillitis tissue. Positive WB detected in HL60 cells, Jurkat cells. Observed molecular weight by Western-blot: 21 kDa

## Preparation

This antibody was obtained by immunization of RAC2 recombinant protein (Accession Number: NM\_002872). Purification method: Protein G purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:107495(RAC2 antibody) at dilution of 1:500



Immunohistochemistry of paraffin-embedded human spleen slide using Catalog No:107495(RAC2 Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human spleen slide using Catalog No:107495(RAC2 Antibody) at dilution of 1:50