

# COX6B2 antibody



Catalog Number: 109495

## Product name

COX6B2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

COX6B2 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue. Positive IHC detected in mouse testis tissue. Observed molecular weight by Western-blot: 11kd

## Preparation

This antibody was obtained by immunization of COX6B2 recombinant protein (Accession Number: BC064548). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


## Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:109495(COX6B2 antibody) at dilution of 1:100



Immunohistochemistry of paraffin-embedded mouse testis tissue slide using Catalog No:109495(COX6B2 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded mouse testis tissue slide using Catalog No:109495(COX6B2 Antibody) at dilution of 1:50 (under 40x lens)