

# Multi tags antibody



Catalog Number: 112895

## Product name

Multi tags antibody

## Specificity

Human; other species not tested.

## Antibody description

Multi tags Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

## Preparation

This antibody was obtained by immunization of Multi tags recombinant protein. Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:5000-1:50000

## Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:112895(V5 tag Antibody) at dilution of 1:10000