

Anti-SIRPB1 antibody

Catalog Number: 175895

Product name

Polyclonal

Anti-SIRPB1 antibody

Ig Type**Specificity**

Rabbit IgG

Human

Applications**Antibody description**

WB

Rabbit polyclonal antibody to SIRPB1

Dilutions**Preparation**

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human sirpb1. aa range:231-280

Validations**Formulation**

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.


Clonality

Fig1.; Western Blot analysis of various cells using CD172b Polyclonal Antibody