

TBXAS1 antibody



Catalog Number: 115895

Product name

TBXAS1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TBXAS1 Rabbit Polyclonal antibody. Positive WB detected in A549 cells, human colon tissue. Positive IHC detected in human lung cancer tissue. Observed molecular weight by Western-blot: 60 kDa

Preparation

This antibody was obtained by immunization of TBXAS1 recombinant protein (Accession Number: NM_001166254). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:115895(TBXAS1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:115895(TBXAS1 antibody) at dilution of 1:50 (under 10x lens)