

# Anti-IKB epsilon Rabbit antibody



Catalog Number: 167895

## Product name

Rabbit IgG

Anti-IKB epsilon Rabbit antibody

## Applications

## Specificity

WB, IHC-P, ICC/IF, IP

Human

## Dilutions

## Antibody description

WB: 1/1000

IKB epsilon Rabbit monoclonal antibody

IHC: 1/50

## Preparation

ICC/IF: 1/50

Antigen: A synthetic peptide of human IKB epsilon

IP: 1/20

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Validations

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Western blot detection of IKB epsilon in HeLa,Ramos using IKB epsilon Rabbit mAb(1:1000 diluted)

