

Anti-Glypican 1/GPC1 Rabbit antibody



Catalog Number: 168895

Product name

Anti-Glypican 1/GPC1 Rabbit antibody

Specificity

Human

Antibody description

Glypican 1/GPC1 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Glypican 1/GPC1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC


Dilutions

1:1000

1:100

1:20

Validations



Western blot detection of Glypican 1/GPC1 in Jurkat, HeLa cell lysates using Glypican 1/GPC1 Rabbit mAb(1:1000 diluted). Predicted band size: 62kDa. Observed band size: 55kDa.