

Human PVR/PVS (His tag) recombinant protein



Catalog Number: 519895

General Information

Gene Name Synonym

Poliovirus Receptor, Nectin-Like Protein 5, NECL-5, CD155, PVR, PVS

Protein Construction

A DNA sequence encoding the Human PVR/PVS Trp21-Asn343 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human PVR/PVS (His tag) recombinant protein