

ARR3 Polyclonal Antibody



Catalog Number: 164195

Product name

ARR3 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to ARR3

Preparation

Antigen: A synthetic peptide of human ARR3

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:200 - 1:2000

Validations



Western blot - ARR3 Polyclonal Antibody

Western blot analysis of extracts of mouse brain, using ARR3 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.