

# ESR2 antibody



Catalog Number: 107195

## Product name

ESR2 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ESR2 Mouse Monoclonal antibody. Positive WB detected in MCF-7 cells, human brain tissue, SW 1990 cells. Observed molecular weight by Western-blot: 50 kDa

## Preparation

This antibody was obtained by immunization of ESR2 recombinant protein (Accession Number: XM\_047431078). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:107195(ESR2 Antibody) at dilution of 1:2000