

Anti-CD19 (3G7) Mouse antibody



Catalog Number: 168195

Product name

Anti-CD19 (3G7) Mouse antibody

Specificity

Human

Antibody description

CD19 (3G7) Mouse monoclonal antibody

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

ICC/IF, FC


Dilutions

ICC: 1/100

Validations



Flow cytometric analysis of Jurkat T cells (black) and Ramos B cells (red), using anti-CD19 mouse mAb.



Confocal immunofluorescent analysis of Ramos (positive cell, left) and Jurkat (negative cell, right) using anti-CD19 mouse mAb (dilution 1:100).